

THE EYE DROP LABEL

Last Updated: Tuesday, 3 February 2009

Author: Karen Stancliffe


Summary	1
Details	1

Summary


Eye drop containers are much smaller than most medication packaging. Often the main label is too big to fit on the bottle. The Fred5 label format allows for a special eye drop label to be printed for use on eye drop medication. The label can also be configured to print for other types of medication.

Details

The setting for label size is found in the menu - Lists > Drugs as shown below.




The label size can be set to S for small or L for Large. By default all eye drop medications are set to S for small. This means an eye drop label will print for any eye drop medication automatically.






To set the label size to S for small for another medication –

- Click on the Lists menu and choose Drugs
- Press F for Find and type in the name of the item you are searching for



Change the Label Size setting from L to S for Small Eye Drop Label.






A label that prints for the small label size setting will print as follows -



A label that prints for the large Label size will print as follows -



PFRED.HELP

1300 731 888 (Local Call Cost) help@fred.com.au

8.30am - 9.00pm (EST) Mon - Fri

8.30am - 5.00pm (EST) Sat, Sun & Public Holidays

www.fred.com.au

A PCA NU Systems Company

