

## Printing long directions on multiple labels explained.

Last Updated: Wednesday, 8 October 2008

Author: Karen Stancliffe

|               |   |
|---------------|---|
| Summary ..... | 1 |
| Details ..... | 1 |

### Summary

The Lockie Thermal (L4) Label only has room for two lines of directions. Even when a compressed font is used, not a lot of text will fit on the label.

If the directions (even with compressed print) will not fit on one label, Fred Dispense recognises that 2 labels must be used. Fred Dispense prints both the labels in the larger (non-compressed) font. This is because even with the compressed print, more than one label is needed.

### Details

For example,

*Take ONE tablet TWICE a day Take on empty stomach, at least half an hour before meals and at bedtime`*

These directions, even when compressed, require 2 labels to print - the last 3 words `and at bedtime` would go onto the 2nd label. Since Fred Dispense has to use 2 labels anyway, the larger non-compressed font is used.

*Take ONE tablet TWICE a day Take on empty stomach, at least half an hour before meals`*

(3 words less than the previous example) can fit in condensed font onto one label - so Fred Dispense would use the condensed font to print this onto one label.

The same functionality is also used for the F5 Thermal Label printing, but the F5 Label has 3 lines of directions, so when compressed this label will fit a lot more.

### FRED.HELP

1300 731 888 (Local Call Cost)

help@fred.com.au

8.30am - 9.00pm (EST) Mon – Fri

8.30am - 5.00pm (EST) Sat, Sun & Public Holidays

[www.fred.com.au](http://www.fred.com.au)

A PCA NU Systems Company